

PKA C gamma Rabbit mAb

Catalog No: #57321

Package Size: #57321-1 50ul #57321-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PKA C gamma Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human PRKACG
Other Names	PRKACG; cAMP-dependent protein kinase catalytic subunit gamma; PKA C-gamma
Accession No.	Swiss-Prot:P22612 Gene ID:5568
Uniprot	P22612
GeneID	5568;
Calculated MW	40 kDa
SDS-PAGE MW	40 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of PRKACG in mouse testis, mouse lung lysates using PKA C gamma antibody.



Western blot analysis of PRKACG in MCF-7, C6, HeLa lysates using PRKACG antibody.

Note: This product is for in vitro research use only