

PRSS8 Rabbit mAb

Catalog No: #57324

Package Size: #57324-1 50ul #57324-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PRSS8 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human PRSS8
Other Names	CAP1; PROSTASIN
Accession No.	Swiss-Prot:Q16651 Gene ID:5652
Uniprot	Q16651
GeneID	5652;
Calculated MW	36 kDa
SDS-PAGE MW	36 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of PRSS8 in MCF-7 lysates using PRSS8 antibody.

Note: This product is for in vitro research use only