

PSF1 Rabbit mAb

Catalog No: #57325

Package Size: #57325-1 50ul #57325-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PSF1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human PSF1
Other Names	GINS complex subunit 1
Accession No.	Swiss-Prot:Q14691 Gene ID:9837
Uniprot	Q14691
GeneID	9837;
Calculated MW	23 kDa
SDS-PAGE MW	23 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

Images



Western blot analysis of PSF1 in rat Brain, C6, 3T3, HeLa lysates using PSF1 antibody.

Immunocytochemistry analysis of PSF1 (green) in Jurkat using PSF1 antibody, and DAPI(blue).



Note: This product is for in vitro research use only