

## PTN Rabbit mAb

Catalog No: #57328

Package Size: #57328-1 50ul #57328-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PTN Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human PTN
Other Names	PTN; HBNF1; NRG1; Pleiotrophin; PTN; Heparin-binding brain mitogen; HBBM; Heparin-binding growth factor 8; HBGF-8; Heparin-binding growth-associated molecule; HB-GAM; Heparin-binding neurite outgrowth-promoting factor 1; HBNF-1; Osteoblast-specific factor 1; OSF-1
Accession No.	Swiss-Prot:P21246 Gene ID:5764
Uniprot	P21246
GeneID	5764;
Calculated MW	19 kDa
SDS-PAGE MW	19 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

IP:1/20WB: 1/500-1/1000

## Images



Western blot analysis of PTN in U87-MG lysates using PTN antibody.



Western blot analysis of PTN in C6, rat Brain lysates using PTN antibody.

---

Note: This product is for in vitro research use only