

## PUS1 Rabbit mAb

Catalog No: #57329



Package Size: #57329-1 50ul #57329-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PUS1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Hamster,Rat
Immunogen Description	Recombinant protein of human PUS1
Other Names	MLASA1
Accession No.	Swiss-Prot:Q9Y606 Gene ID:80324
Uniprot	Q9Y606
GeneID	80324;
Calculated MW	47 kDa
SDS-PAGE MW	44,47 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

Note: This product is for in vitro research use only