

RALB Rabbit mAb

Catalog No: #57331

Package Size: #57331-1 50ul #57331-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RALB Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human RALB
Other Names	Ras-related protein Ral-B
Accession No.	Swiss-Prot:P11234 Gene ID:5899
Uniprot	P11234
GeneID	5899;
Calculated MW	23 kDa
SDS-PAGE MW	23 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of RALB in mouse colon lysates using RALB antibody.

Western blot analysis of RALB in MCF-7 lysates using RALB antibody.



Note: This product is for in vitro research use only