

S100 alpha9 Rabbit mAb

Catalog No: #57334

Package Size: #57334-1 50ul #57334-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	S100 alpha9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of mouse S100A9
Other Names	S100 alpha9; Calgranulin-B; Leukocyte L1 complex heavy chain
Accession No.	Swiss-Prot:P31725 Gene ID:20202
Uniprot	P31725
GeneID	20202;
Calculated MW	13 kDa
SDS-PAGE MW	13 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200

Images



Western blot analysis of S100A9 in mouse spleen, mouse lung lysates using S100 alpha9 antibody.

Note: This product is for in vitro research use only