

## SODD Rabbit mAb

Catalog No: #57340

Package Size: #57340-1 50ul #57340-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SODD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human Silencer of Death Domain
Other Names	BAG4; SODD; BAG family molecular chaperone regulator 4; BAG-4; Bcl-2-associated athanogene 4; Silencer of death domains
Accession No.	Swiss-Prot:O95429 Gene ID:9530
Uniprot	O95429
GeneID	9530;
Calculated MW	50 kDa
SDS-PAGE MW	70 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

## Images



Western blot analysis of Silencer of Death Domain in MOLT4 lysates using SODD antibody.

Note: This product is for in vitro research use only