

SLBP Rabbit mAb

Catalog No: #57342

Package Size: #57342-1 50ul #57342-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human SLBP
Other Names	HBP
Accession No.	Swiss-Prot:Q14493 Gene ID:7884
Uniprot	Q14493
GeneID	7884;
Calculated MW	31 kDa
SDS-PAGE MW	31 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of SLBP in 293T, Jurkat, Ramos lysates using SLBP antibody.

Note: This product is for in vitro research use only