

SLC25A12 Rabbit mAb

Catalog No: #57343

Package Size: #57343-1 50ul #57343-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SLC25A12 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Mouse,Rat
Immunogen Description	A synthetic peptide of human SLC25A12
Other Names	AGC1; ARALAR; EIEE39
Accession No.	Swiss-Prot:O75746 Gene ID:8604
Uniprot	O75746
GeneID	8604;
Calculated MW	75 kDa
SDS-PAGE MW	75 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of SLC25A12 in rat Brain, C6, 3T3 lysates using SLC25A12 antibody.

Note: This product is for in vitro research use only