

Striatin 4 Rabbit mAb

Catalog No: #57347

Package Size: #57347-1 50ul #57347-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Striatin 4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human Striatin 4
Other Names	ZIN; PPP2R6C; zinedin
Accession No.	Swiss-Prot:Q9NRL3 Gene ID:29888
Uniprot	Q9NRL3
GeneID	29888;
Calculated MW	81 kDa
SDS-PAGE MW	81 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000ICC: 1/50-1/200

Images



Western blot analysis of Striatin 4 in HeLa lysates using Striatin 4 antibody.

Note: This product is for in vitro research use only