

TACSTD2 Rabbit mAb

Catalog No: #57348

Package Size: #57348-1 50ul #57348-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TACSTD2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human TACD2
Other Names	TACSTD2; GA733-1; M1S1; TROP2; Tumor-associated calcium signal transducer 2; Cell surface glycoprotein Trop-2; Membrane component chromosome 1 surface marker 1; Pancreatic carcinoma marker protein GA733-1
Accession No.	Swiss-Prot:P09758 Gene ID:4070
Uniprot	P09758
GeneID	4070;
Calculated MW	36 kDa
SDS-PAGE MW	45-65 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of TROP2 in HCT116, PC-3, MCF-7 lysates using TACSTD2 antibody.

Note: This product is for in vitro research use only