

TAF1C Rabbit mAb

Catalog No: #57349

Package Size: #57349-1 50ul #57349-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TAF1C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human TAF1C
Other Names	SL1; TAFI95; TAFI110; MGC:39976
Accession No.	Swiss-Prot:Q15572 Gene ID:9013
Uniprot	Q15572
GeneID	9013;
Calculated MW	95 kDa
SDS-PAGE MW	95 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of TAF1C in K562 lysates using TAF1C antibody.

Note: This product is for in vitro research use only