

## TIAM2 Rabbit mAb

Catalog No: #57350

Package Size: #57350-1 50ul #57350-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TIAM2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human TIAM2
Other Names	STEF; TIAM-2
Accession No.	Swiss-Prot:Q8IVF5 Gene ID:26230
Uniprot	Q8IVF5
GeneID	26230;
Calculated MW	190 kDa
SDS-PAGE MW	190 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200

## Images



Western blot analysis of TIAM2 in K562, 293 lysates using TIAM2 antibody.

Note: This product is for in vitro research use only