

## TIMM44 Rabbit mAb

Catalog No: #57351

Package Size: #57351-1 50ul #57351-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TIMM44 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human TIMM44
Other Names	TIM44
Accession No.	Swiss-Prot:O43615 Gene ID:10469
Uniprot	O43615
GeneID	10469;
Calculated MW	51 kDa
SDS-PAGE MW	45 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details


IP: 1/20WB: 1/500-1/1000ICC: 1/50-1/200

## Images



Western blot analysis of TIMM44 in HeLa lysates using TIMM44 antibody.

Immunocytochemistry analysis of TIMM44(green) in Hela using TIMM44 antibody, and DAPI(blue)



---

Note: This product is for in vitro research use only