

## Transmembrane Protein 161B Rabbit mAb

Catalog No: #57352

Package Size: #57352-1 50ul #57352-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Transmembrane Protein 161B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human TMEM161B
Other Names	Transmembrane protein 161B
Accession No.	Swiss-Prot:Q8NDZ6 Gene ID:153396
Uniprot	Q8NDZ6
GeneID	153396;
Calculated MW	55 kDa
SDS-PAGE MW	55 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP:1/20WB: 1/500-1/1000

## Images



Western blot analysis of TMEM161B in K562, HeLa lysates using Transmembrane Protein 161B antibody.

Note: This product is for in vitro research use only