

Thrombopoietin Receptor Rabbit mAb

Catalog No: #57353

Package Size: #57353-1 50ul #57353-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Thrombopoietin Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human TPOR
Other Names	MPL; TPOR; Thrombopoietin receptor; TPO-R; Myeloproliferative leukemia protein; Proto-oncogene c-Mpl; CD110
Accession No.	Swiss-Prot:P40238 Gene ID:4352
Uniprot	P40238
GeneID	4352;
Calculated MW	71 kDa
SDS-PAGE MW	85 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000

Images



Western blot analysis of TPOR/MPL in K562, C6 lysates using Thrombopoietin Receptor antibody.

Note: This product is for in vitro research use only