

TRIF Rabbit mAb

Catalog No: #57354

Package Size: #57354-1 50ul #57354-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRIF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human TRIF
Other Names	TRIF; IIAE6; MyD88-3; PRVTIRB; TICAM-1
Accession No.	Swiss-Prot:Q8IUC6 Gene ID:148022
Uniprot	Q8IUC6
GeneID	148022;
Calculated MW	76 kDa
SDS-PAGE MW	98 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of TRIF in rat testis lysates using TRIF antibody.



Western blot analysis of TRIF in K562, 3T3, HeLa lysates using TRIF antibody.

Note: This product is for in vitro research use only