

TRIM21 Rabbit mAb

Catalog No: #57355

Package Size: #57355-1 50ul #57355-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TRIM21 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human TRIM21
Other Names	SSA; RO52; SSA1; RNF81; Ro/SSA
Accession No.	Swiss-Prot:P19474 Gene ID:6737
Uniprot	P19474
GeneID	6737;
Calculated MW	54 kDa
SDS-PAGE MW	50 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000

Images



Western blot analysis of TRIM21 in K562, HeLa lysates using TRIM21 antibody.

Note: This product is for in vitro research use only