

## Transcription Termination Factor 2 Rabbit mAb

Catalog No: #57356

Package Size: #57356-1 50ul #57356-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Transcription Termination Factor 2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human TTF2
Other Names	HuF2; ZGRF6
Accession No.	Swiss-Prot:Q9UNY4 Gene ID:8458
Uniprot	Q9UNY4
GeneID	8458;
Calculated MW	130 kDa
SDS-PAGE MW	130 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200


## Images



Western blot analysis of TTF2 in K562, rat Brain, C6 lysates using Transcription Termination Factor 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TTF2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TTF2 (green) in HeLa using TTF2 antibody, and DAPI (blue)

---

Note: This product is for in vitro research use only