

U1C Rabbit mAb

Catalog No: #57358

Package Size: #57358-1 50ul #57358-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	U1C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human U1-C
Other Names	U1 small nuclear ribonucleoprotein C; U1 snRNP C; U1-C
Accession No.	Swiss-Prot:P09234 Gene ID:6631
Uniprot	P09234
GeneID	6631;
Calculated MW	17 kDa
SDS-PAGE MW	21 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of U1C in rat lung lysates using U1C antibody.



Western blot analysis of U1C in THP-1, Jurkat, HeLa, K562 lysates using U1C antibody.

Note: This product is for in vitro research use only