

UCHL5 Rabbit mAb

Catalog No: #57360

Package Size: #57360-1 50ul #57360-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	UCHL5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human UCH37
Other Names	UCH37; CGI-70; INO80R; UCH-L5
Accession No.	Swiss-Prot:Q9Y5K5 Gene ID:51377
Uniprot	Q9Y5K5
GeneID	51377;
Calculated MW	38 kDa
SDS-PAGE MW	38 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP: 1/20WB: 1/500-1/1000ICC: 1/50-1/200


Images



Western blot analysis of UCH37 in rat spleen lysates using UCHL5 antibody.



Western blot analysis of UCH37 in 3T3, HeLa lysates using UCH37 antibody.



Immunocytochemistry analysis of UCH37(green) in HeLa using UCH37 antibody, and DAPI(blue)

Note: This product is for in vitro research use only