

## Complex III Subunit 7 Rabbit mAb

Catalog No: #57361

Package Size: #57361-1 50ul #57361-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Complex III Subunit 7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of human UQCRB
Other Names	QPC; QCR7; QP-C; UQBC; UQBP; UQPC; UQCR6; MC3DN3
Accession No.	Swiss-Prot:P14927 Gene ID:7381
Uniprot	P14927
GeneID	7381;
Calculated MW	14 kDa
SDS-PAGE MW	14 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

## Images



Western blot analysis of UQCRB in rat Brain lysates using Complex III Subunit 7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UQCRB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

---

Note: This product is for in vitro research use only