

ZNF160 Rabbit mAb

Catalog No: #57363

Package Size: #57363-1 50ul #57363-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF160 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ZNF160
Other Names	Zinc finger protein 160; Zinc finger protein HZF5; Zinc finger protein Kr18
Accession No.	Swiss-Prot:Q9HCG1 Gene ID:90338
Uniprot	Q9HCG1
GeneID	90338;
Calculated MW	94 kDa
SDS-PAGE MW	74 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of ZNF160 in mouse liver lysates using ZNF160 antibody.

Note: This product is for in vitro research use only