

FIG4 Rabbit mAb

Catalog No: #57366

Package Size: #57366-1 50ul #57366-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FIG4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human FIG4
Other Names	YVS; BTOP; SAC3; ALS11; CMT4J; KIAA0274; dJ249I4.1
Accession No.	Swiss-Prot:Q92562 Gene ID:9896
Uniprot	Q92562
GeneID	9896;
Calculated MW	104 kDa
SDS-PAGE MW	104 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of FIG4 in HL-60, U2OS, U87-MG lysates using FIG4 antibody.

Note: This product is for in vitro research use only