

Activin Receptor Type IIB Rabbit mAb

Catalog No: #57368

Package Size: #57368-1 50ul #57368-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Activin Receptor Type IIB Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB
Species Reactivity	Human,Mouse,Rat,Hamster
Immunogen Description	A synthetic peptide of human Activin Receptor Type IIB/ACVR2B
Other Names	ActR IIB; ACVR 2B; HTX4
Accession No.	Swiss-Prot:Q13705 Gene ID:93
Uniprot	Q13705
GeneID	93;
Calculated MW	58 kDa
SDS-PAGE MW	58 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000

Images



Western blot analysis of Activin Receptor Type IIB/ACVR2B in CHO-K1, mouse Brain lysates using Activin Receptor Type IIB antibody.

Note: This product is for in vitro research use only