

ACTR1B Rabbit mAb

Catalog No: #57369

Package Size: #57369-1 50ul #57369-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACTR1B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ACTR1B
Other Names	PC3; ARP1B; CTRN2
Accession No.	Swiss-Prot:P42025 Gene ID:10120
Uniprot	P42025
GeneID	10120;
Calculated MW	42 kDa
SDS-PAGE MW	42 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200


Images



Western blot analysis of ACTR1B in HeLa, A549, HL-60, U251, U87-MG lysates using ACTR1B antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ACTR1B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of ACTR1B (green) in HeLa using ACTR1B antibody, and DAPI (blue)

Note: This product is for in vitro research use only