

## Actin Related Protein 2 Rabbit mAb

Catalog No: #57370

Package Size: #57370-1 50ul #57370-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Actin Related Protein 2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Arp2
Other Names	ACTR2; ARP2; ARP2/3 COMPLEX
Accession No.	Swiss-Prot:P61160 Gene ID:10097
Uniprot	P61160
GeneID	10097;
Calculated MW	45 kDa
SDS-PAGE MW	45 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/500-1/1000 ICC: 1/50-1/200

## Images



Western blot analysis of Arp2 in C6, mouse Brain, Ramos lysates using Arp2 antibody.



Immunocytochemistry analysis of Actin Related Protein 2 (green) in HeLa using Actin Related Protein 2 antibody, and DAPI(blue).

---

Note: This product is for in vitro research use only