

FLT3L Antibody FITC Conjugated

Catalog No: #C00634F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FLT3L Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-120 235 derived from human FLT3L
Conjugates	FITC
Target Name	FLT3L
Other Names	FL; FLT3L; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; SL cytokine; FLT3LG
Accession No.	Swiss-Prot#P49771NCBI Gene ID2323
Uniprot	P49771
GeneID	2323;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Stimulates the proliferation of early hematopoietic cells by activating FLT3. Synergizes well with a number of other colony stimulating factors and interleukins.

Note: This product is for in vitro research use only