

GMEB1 Rabbit mAb

Catalog No: #57376

Package Size: #57376-1 50ul #57376-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GMEB1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human GMEB1
Other Names	PIF96; P96PIF
Accession No.	Swiss-Prot:Q9Y692 Gene ID:10691
Uniprot	Q9Y692
GeneID	10691;
Calculated MW	63 kDa
SDS-PAGE MW	80 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of GMEB1 in A549, HL-60, U2OS lysates using GMEB1 antibody.

Note: This product is for in vitro research use only