

IL-1 Receptor Accessory Protein Rabbit mAb

Catalog No: #57378

Package Size: #57378-1 50ul #57378-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL-1 Receptor Accessory Protein Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human IL-1RAcP
Other Names	IL-1 Receptor Accessory Protein; IL1R3; C3orf13; IL-1RAcP
Accession No.	Swiss-Prot:Q9NPH3 Gene ID:3556
Uniprot	Q9NPH3
GeneID	3556;
Calculated MW	65 kDa
SDS-PAGE MW	78 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of IL-1RAcP in HL-60, U251, U87-MG lysates using IL-1 Receptor Accessory Protein antibody.

Note: This product is for in vitro research use only