

## LILRB1 Rabbit mAb

Catalog No: #57380

Package Size: #57380-1 50ul #57380-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	LILRB1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human LILRB1
Other Names	ILT2; LIR1; MIR7; PIRB; CD85J; ILT-2; LIR-1; MIR-7; PIR-B
Accession No.	Swiss-Prot:Q8NHL6 Gene ID:10859
Uniprot	Q8NHL6
GeneID	10859;
Calculated MW	71 kDa
SDS-PAGE MW	110 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP:1/20WB: 1/500-1/1000

## Images



Western blot analysis of LILRB1 in HL-60, U251 lysates using LILRB1 antibody.

Note: This product is for in vitro research use only