

MYH9 Rabbit mAb

Catalog No: #57384

Package Size: #57384-1 50ul #57384-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MYH9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human non-muscle Myosin IIA
Other Names	MHA; FTNS; EPSTS; BDPLT6; DFNA17; MATINS; NMMHCA; NMHC-II-A; NMMHC-IIA
Accession No.	Swiss-Prot:P35579 Gene ID:4627
Uniprot	P35579
GeneID	4627;
Calculated MW	227 kDa
SDS-PAGE MW	227 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of nonmuscle Myosin IIA in HeLa, A549, HL-60, U2OS, C6 lysates using MYH9 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using nonmuscle Myosin IIA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note: This product is for in vitro research use only