

PTP4A2 Rabbit mAb

Catalog No: #57389

Package Size: #57389-1 50ul #57389-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PTP4A2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human PTP4A2/PRL2
Other Names	protein tyrosine phosphatase type IVA; member 2; HH13; OV-1; PRL2; HH7-2; PRL-2; PTP4A; HU-PP-1; PTPCAAX2; ptp-IV1a; ptp-IV1b
Accession No.	Swiss-Prot:Q12974 Gene ID:8073
Uniprot	Q12974
GeneID	8073;
Calculated MW	19 kDa
SDS-PAGE MW	19 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of PTP4A2/PRL2 in HeLa, A549, HL-60 lysates using PTP4A2 antibody.

Note: This product is for in vitro research use only