

QSOX1 Rabbit mAb

Catalog No: #57390

Package Size: #57390-1 50ul #57390-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	QSOX1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human Quiescin Q6
Other Names	Q6; QSCN6
Accession No.	Swiss-Prot:O00391 Gene ID:5768
Uniprot	O00391
GeneID	5768;
Calculated MW	83 kDa
SDS-PAGE MW	83 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Quiescin Q6 in A549, U2OS lysates using QSOX1 antibody.

Note: This product is for in vitro research use only