

Semaphorin 7A Rabbit mAb

Catalog No: #57392

Package Size: #57392-1 50ul #57392-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Semaphorin 7A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human Semaphorin 7a
Other Names	JMH; CD108; SEMAL; CDw108; SEMAK1; H-Sema-L; H-SEMA-K1
Accession No.	Swiss-Prot:O75326 Gene ID:8482
Uniprot	O75326
GeneID	8482;
Calculated MW	75 kDa
SDS-PAGE MW	75 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of Semaphorin 7a in A549, HL-60, U251, U87-MG lysates using Semaphorin 7A antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Semaphorin 7a antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note: This product is for in vitro research use only