

## WVOX Rabbit mAb

Catalog No: #57398

Package Size: #57398-1 50ul #57398-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	WVOX Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human WVOX
Other Names	FOR; WOX1; EIEE28; FRA16D; SCAR12; HHCMA56; PRO0128; SDR41C1; D16S432E
Accession No.	Swiss-Prot:Q9NZC7 Gene ID:51741
Uniprot	Q9NZC7
GeneID	51741;
Calculated MW	47 kDa
SDS-PAGE MW	47 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100


## Images



Western blot analysis of WVOX in mouse testis, mouse spleen lysates using WVOX antibody.



Western blot analysis of WWOX in HeLa, U2OS, U251MG lysates using WWOX antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using WWOX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

---

Note: This product is for in vitro research use only