

YTHDC2 Rabbit mAb

Catalog No: #57399

Package Size: #57399-1 50ul #57399-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	YTHDC2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of mouse YTHDC2
Other Names	CAHL; hYTHDC2
Accession No.	Swiss-Prot:B2RR83 Gene ID:240255
Uniprot	B2RR83
GeneID	240255;
Calculated MW	161 kDa
SDS-PAGE MW	161 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of YTHDC2 in HeLa, A549, HL-60, U251, U87-MG lysates using YTHDC2 antibody.

Note: This product is for in vitro research use only