

# Mouse Low Density Lipoprotein Receptor Related Protein 6 (LRP6) ELISA Kit

Catalog No: #EK20580

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20580-1 48T #EK20580-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Low Density Lipoprotein Receptor Related Protein 6 (LRP6) ELISA Kit
Brief Description	ELISA Kit
Host Species	Mouse
Purification	ELISA

## Application Details

Detect Range:0.156-10ng/ml

Sensitivity:0.067ng/mL


Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Images



## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only