

## ALKBH5 Polyclonal Antibody

Catalog No: #48826

Package Size: #48826-1 50ul #48826-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description


Product Name	ALKBH5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of human ALKBH5
Immunogen Description	Synthesized peptide derived from human ALKBH5
Target Name	ALKBH5 ABH5 OFOX1
Other Names	Probable alpha-ketoglutarate-dependent dioxygenase ABH5 (EC 1.14.11.-.;Alkylated DNA repair protein alkB homolog 5)
Uniprot	Q6P6C2
GeneID	54890
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

IHC-p 1:50-200, WB 1:500-2000. IF 1:50-200

## Images





Immunohistochemical analysis of paraffin-embedded human tonsil.

---

Note: This product is for in vitro research use only