

TRPV3 Polyclonal Antibody

Catalog No: #48836

Package Size: #48836-1 50ul #48836-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRPV3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	TRPV3 Polyclonal Antibody detects endogenous levels of TRPV3 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human TRPV3.
Target Name	TRPV3
Other Names	TRPV3; Transient receptor potential cation channel subfamily V member 3; TrpV3; Vanilloid receptor-like 3; VRL-3
Uniprot	Q8NET8
GeneID	162514
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC-p: 1:100-300 ELISA: 1:20000.. IF 1:50-200

Images



Western Blot analysis of extracts from CACO-2 cells, using TRPV3 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:100

Note: This product is for in vitro research use only