

Human Low Density Lipoprotein Receptor (LDLR) ELISA Kit

Catalog No: #EK20583

Orders: order@signalwayantibody.com

Package Size: #EK20583-1 48T #EK20583-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Low Density Lipoprotein Receptor (LDLR) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.118ng/mL


Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Images



Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only