

ASGL1 rabbit pAb

Catalog No: #48840

Package Size: #48840-1 50ul #48840-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ASGL1 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of ASGL1 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human ASGL1
Target Name	ASRGL1 ALP CRASH
Uniprot	Q7L266
GeneID	80150
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1 μ g 500-2000

Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only