

Rat NACHT, LRR and PYD Domains-Containing Protein 3 (NLRP3) ELISA Kit

Catalog No: #EK20596

Orders: order@signalwayantibody.com

Package Size: #EK20596-1 48T #EK20596-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Rat NACHT, LRR and PYD Domains-Containing Protein 3 (NLRP3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.054ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only