

Human Nuclear receptor subfamily 1 group D member 1 (NR1D1) ELISA Kit



Catalog No: #EK20609

Orders: order@signalwayantibody.com

Package Size: #EK20609-1 48T #EK20609-2 96T

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human Nuclear receptor subfamily 1 group D member 1 (NR1D1) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Human |

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:24.5pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only