

Recombinant Human Platelet-derived Growth Factor-AA



Catalog No: #AP60081

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP60081-1 10ug #AP60081-2 100ug #AP60081-3 500ug

Description

Product Name	Recombinant Human Platelet-derived Growth Factor-AA
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE.
Calculated MW	Approximately 28.8 kDa, a disulfide-linked homodimeric protein containing two 126 amino acid residues polypeptide. But it migrates with an apparent molecular mass of 33.6 kDa in SDS-PAGE.
Target Sequence	MSIEEAVPAV CKTRTVIYEI PRSQVDPTSA NFLIWPPCVE VKRCTGCCNT SSVKQCPSRV HHRSVKVAKV EYVRKKPKLK EVQVRLEEHL ECACATTSLN PDYREEDTGR PRESGKKRKR KRLKPT
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB, 300 mM NaCl, pH 6.0.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only