

Recombinant Human Placenta Growth Factor-2

Catalog No: #AP60090



Package Size: #AP60090-1 5ug #AP60090-2 100ug #AP60090-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Placenta Growth Factor-2
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	P49763-3
GeneID	5228
Calculated MW	Approximately 34.6 kDa, a disulfide-linked homodimer consisting of two 152 amino acid polypeptide chains.
Target Sequence	LPAVPPQQWA LSAGNGSSEV EVVPFQEVWG RSYCRALERL VDVVSEYPSE VEHMFSPSCV SLLRCTGCCG DENLHCVPE TANVTMQLLK IRSGDRPSYV ELTFSQHVRC ECRPLREKMK PERRRPKGRG KRRREKQRPT DCHLCGDAVP RR
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, with 0.02 % Tween-20.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only