

Recombinant Human Neuregulin 1-beta1 EGF Domain, 177-241a.a.



Catalog No: #AP60109

Orders: order@signalwayantibody.com

Package Size: #AP60109-1 10ug #AP60109-2 100ug #AP60109-3 500ug

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Neuregulin 1-beta1 EGF Domain, 177-241a.a.
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	Q02297
GeneID	3084
Calculated MW	Approximately 7.5 kDa, a single non-glycosylated polypeptide chain containing 65 amino acids.
Target Sequence	SHLVKCAEKE KTFCVNGGEC FMVKDLSNPS RYLCKPNEF TGDRQCQNYVM ASFYKHLGIE FMEAE
Formulation	Lyophilized from a 0.2 µm filtered solution in 1 x PBS, pH 7.4, with 5 % trehalose.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only