

Recombinant Human Neurturin

Catalog No: #AP60119



Package Size: #AP60119-1 5ug #AP60119-2 100ug #AP60119-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Neurturin
Host Species	Escherichia coli.
Purification	> 96 % by SDS-PAGE and HPLC analyses.
Uniprot	Q99748
GeneID	4902
Calculated MW	Approximately 23.4 kDa, a homodimeric protein consisting of two 102 amino acid non-glycosylated polypeptide chains.
Target Sequence	ARLGARPCGL RELEVRVSEL GLGYASDETV LFRYCAGACE AAARVYDLGL RRLRQRRRLR RERVRAQPCC RPTAYEDEVSLFLDAHSRYHT VHELARECA CV
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 30 mM Citrate Sodium, pH 4.2, 400 mM NaCl, with 0.02 % Tween-20.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only