

# Recombinant Rhesus Macaque Granulocyte Colony Stimulating Factor

Catalog No: #AP60146


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP60146-1 10ug #AP60146-2 100ug #AP60146-3 500ug

## Description

Product Name	Recombinant Rhesus Macaque Granulocyte Colony Stimulating Factor
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Uniprot	F7H1Q6
GeneID	698961
Calculated MW	Approximately 18.9 kDa, a single non-glycosylated polypeptide chain containing 174 amino acids.
Target Sequence	TPLGPASSLP QSFLKCLEQ VRKIQGDGAA LQEKLCA <sup>TYK</sup> LCHPEELVLL RHSLGIPWAP LSSCPSQALQ LTGCLSQLHS SLFLYQGLLQ ALEGISPELS PTLDTLQLDI ADFATTIWQQ MEDLGMAPAL QPTQGAMPAF TSAFQRRAGG VLVASHLQRF LELAYRVLRH LAQS
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



Note: This product is for in vitro research use only